



Introducing Abbadox Physician Imaging Portal

Imaging Partners has invested in this upgrade for physicians to access reports and view images more efficiently. Abbadox is a true zero download web viewer available in any browser.



How to access Abbadox Portal

Login using username and password

	bbaDox	
CareElov	w Login	
Username	Login	
Password		
	Log in	
	Forgot Your Password?	
	OR	
	Log in with Single Sign-On (SSO)	



PHYSICIAN RESOURCES

A partnership with Moran, Rowen & Dorsey and St. Joseph Hospital of Orange



Search for patient's images

Type your patient's first and last name.

Be sure the "study performed" field is set to the proper parameters. Then, press enter or click search.







3.

View patient's images

Single-click the "DICOM" Action Button to view patient images

DICOM Button

AbbaDox Careflow	× +					– 0 X
← C ⋒ 🗅 https	://platform.abbadox.com/abba				ame=Worklist%20A	Action;title=;cssClass=null:actionlconBg=%233633E3
						C) Other favorites
III AbbaDox 🗉						> 📮 💿 📅 Test Test
Search	MRN	Patient Last Nametest	Patient First Name	Date of Birth MM/DD/YYYY	Date of Service	Ke Case Id / Accession # Add Filters Search Clear Save
Workflows	Q Search Results					
ttt keterrais >	Date of Service ~	Record Number Patient Name	Referring Physician	Document Visit Des	Actions	
	Showing 1 - 1 of 1	810205 TEST JULIO	NAME PHYSICIAN	Legacy - MG DIG DIAC	e: 75 ~	IMAGING PARTNERS
						PHYSICIAN NAME
						PATIENT: JULIO TEST DOB: 05/01/2012 ID #: 810205 Accession#: 20246/HUO_2QF Facility: Legacy Exam Date: 10/22/2024
					_	PROCEDURE: DIGITAL DIAGNOSTIC 2D MAMMOGRAM BILATERAL WITH CAD AND 3D TOMOSYNTHESIS
						COMPARISON: None.
						INDICATION: Testing
						BREAST COMPOSITION: The breasts are almost entirely fatty.
						FINDINGS: Testing ASSESSMENT
						BI-RADS Category 1 - Negative examination.
						CONCLUSION AND RECOMMENDATIONS: Testing Recommend clinical evaluation.
						Patient was informed of the results.
						PLEASE NOTE: A NORMAL MAMMOGRAM DOES NOT EXCLUDE THE POSSIBILITY OF BREAST CANCER. A CLINICALLY SUSPICIOUS PALPABLE LUMP SHOULD BE BIOPSIED.
						Patient's information is entered into a reminder system with a target due date for the next mammogram.







How to view series of images

To scroll through the series of images you can use your mouse wheel or left click and drag the image number slide bar.

		Image number slide bar
M IPOC-ICA × -	F	
← → C 🌲 pacs.imagingpartner	s.com/webaccess/De	fault.ashx
IBM iConnect Access	ST, JULIO 02-Aug-1990 - 30 yrs old I	Female ID 718016
EXAM LIST HISTORY	×	18 Jan-2019 10:00:00 AM AB DOMEN & PELVIS (STONE PROTOCOL)
One or more studies may be missing or incomplete. Modality Non-DICOI All All	M Content	CT 63 CT 63 CT 62 CT 208 CT 151 CT 61 CT 1 CT 52- 1 100x 53 23x 55 95 55 55 56 58 55 56 56 56 56 56 56 56 56 56 56 56 56
Description Accession#	¥	TEST, JULIO 718/16 F 02 Aug-90 [PR: MARY] 18- an-19 10:00:00 AM
5 out of 5 S	Sort By Newest 👻	1 15 cm
18-Jan-2019 10:00:00 AM CT,Text ABDOMEN & PELVIS (STONE PROTOCOL) Ace: 201926XLFUH7Q		
09-Jul-2015 5:03:00 PM MG DIG SCR BILAT/3D/CAD Acc: 806683		61 CHEST WO SAG CHEST WO SAG
25-Jul-2014 12:21:00 PM MG MAM DIAGNOSTIC DIGITAL WCAD BI Ace: 201452CWJD4A0		FOV: 267 x 316/04 Profile Pro
15-Feb-2012 8:25:00 AM CT CHEST W & ABD WOW & PELVIS W Acc: 201214LEP31JB		718016 A C 718016 F 02-Aug-90 F 02-Aug-90 [PRIMARY] 18-Jan-19 10:00:00 AM 1
02-Aug-2009 8:29:00 AM		







How to use image tools

To access different tools, simply right click anywhere on the image and select your tool preference.





6. I

How to change viewer layout

To alter the image viewer layout, click the layout button and select your preference.







View multiple series of images

To view a different series of images, left click and drag the image series thumbnail to the desired view box.





Turn on/off localizer

You can toggle the appearance of localizer lines on and off by clicking the localizer lines button







9.

How to link image

Some image series can be linked. To link an image series, click the link scroll button and continue to scroll through an image series.



10.

Linked images

The green circle will indicate which images are linked.









View exam report

Click the exam report button to review the radiologist's findings.





How to print report

To print the report, click the print button.







13.

Exit exam and log out

To ensure the privacy of patients, please be sure to exit the exam by closing the browser or sign-Out when finish viewing.



